Spectroscopy Performance Note

Quantitative Depth Profile (QDP) Analysis of Aluminized Steel

- Aluminum Coating Weight
- Coating Thickness
- AI-Fe Alloy Layer Characteristics

Introduction

Aluminum-coated steel products have the benefits of enhanced corrosion resistance, bright metallic appearance, high heat reflectivity and electrical conductivity while retaining the formability of uncoated steel. Additionally the iron-aluminum intermetallic alloy layer provides a second line of defense for oxidation resistance, abrasion resistance and increased hardness.



GDS850A

A specimen of hot dipped aluminized steel was analyzed on the GDS-850A for characterization of the coating. Analysis parameters are shown in Table 1.

Table 1: Method Parameters for Aluminized Steel

Anode Diameter:	4 mm
Lamp Type:	DC
GDS operating Conditions:	40mA, 1000V
Discharge Stabilization:	Current (Control Mode) / Voltage (Pressure Control)
Minimum Data Acquisition Rate:	10 s ⁻¹
Profile Duration:	720 s (thickness dependant)
Cooling:	Cooling, 15°C

An overlay of three replicate analyses is shown in Figure 1. The plot is displayed in concentration in weight percent versus depth in micrometers (μ m). Characteristics of the coating were calculated using GDS-850A software and displayed in Table 2. The table includes the depth of the free aluminum in micrometers (μ m), the thickness of the AI-Fe alloy layer in micrometers (μ m), the coating weight of aluminum in grams per square meter (g/m²) and the concentration of Fe in the coating. Also, the coating weight of C on the surface of the steel is calculated and displayed in grams per square meter (g/m²).




Figure 1. Overlay of three replicate analyses.

Table 2. Coating Characteristics

Name	Free Al, µm	Alloy, µm	Al Coat Wt g/m ²	Fe% Coat	C Wt, interface g/m ²
Aluminized	41.6	12.3	157	2.20	0.038
Aluminized	39.4	12.4	153	2.19	0.037
Aluminized	38.0	12.8	149	2.16	0.037
Mean	39.7	12.5	153	2.18	0.037
Rsd	4.58	2.12	2.61	0.95	1.5

The GDS-850A software can also display variation with depth of selected analytes. The gradients of Al, Fe and C through the sample are shown in Table 3; composition is expressed in weight percent versus depth in micrometers at a granularity of 2 μ m.

For confirmation of the coating thickness the specimen was mounted in cross section and polished. The metallographic cross-sectional image shown at 500X is displayed in Figure 2. The free aluminum and alloy layers of the coating are clearly seen.




Figure 2. Free aluminum and alloy layers of the coating, 500X magnification.

Summary

Quantitative Depth Profile analysis provides coating thickness information comparable with that found by cross-sectional metallography. Simultaneously, the GDS-850A software calculates coating weights, alloy compositions and elemental gradients. Additionally, contamination on the surface or at the interface can be identified and quantified.

Table 3. Gradients of Al, Fe, and C

AI,	Al, Fe, and C						
Depth (µm)		Fe (%)	C (%)				
2	95.7	2.1	0.016				
4	97.3	1.9	0.005				
6	97.5	2.0	0.003				
8	97.5	2.1	0.001				
10	97.5	2.2	0.001				
12	97.6	2.2	0.001				
14	97.5	2.2	0.001				
16	97.6	2.2	0.000				
18	97.6	2.2	0.000				
20	97.6	2.2	0.001				
22	97.7	2.2	0.000				
24	97.6	2.2	0.000				
26	97.5	2.3	0.000				
28	97.3	2.5	0.000				
30	97.0	2.9	0.001				
32	96.3	3.6	0.001				
34	94.2	5.7	0.004				
36	89.7	10.1	0.008				
38	85.0	14.8	0.015				
40	80.5	19.3	0.021				
42	75.9	23.8	0.021				
44	71.0	28.7	0.032				
46	66.0	33.7	0.038				
48	60.6	39.1	0.042				
50	54.7	45.0	0.046				
50	48.9	50.8	0.046				
54	42.6	57.0	0.045				
56	36.7	62.9	0.043				
58	30.8	68.7	0.040				
60	24.8	74.7	0.036				
62	18.3	81.2	0.030				
64	12.1	87.4	0.025				
66	6.5	93.0	0.019				
68	2.4	97.0	0.015				
70	0.7	98.7	0.014				
70	0.2	99.2	0.012				
72	0.2	99.3	0.011				
76		99.4	0.012				
78	0.1	99.4	0.013				
80	0.1	99.4	0.014				
82	0.1	99.3	0.015				
84		99.3					
	0.1	99.3	0.017				
86	0.1	99.3	0.018				
88	0.1						
90 92	0.1	99.3	0.020				
		99.3	0.021				
94	0.1	99.3	0.022				
96	0.1	99.3	0.023				
98	0.1	99.3	0.024				
100	0.1	99.3	0.025				



[®] LECO Corporation
 3000 Lakeview Ave. • St. Joseph, MI 49085
 Phone: 800-292-6141 • Fax: 269-982-8977
 info@leco.com • www.leco.com • ISO-9001:2000 No. FM 24045