


Quantitative Depth Profile (QDP) Analysis of Galvannealed Steel

GDS Applications Lab • LECO Corporation, Saint Joseph, Michigan USA



Name	Depth, μm	Coating wt, g/m^2	Fe%	Al%
Galvanneal	8.06	57.27	9.93	0.19

GDS-850A analysis of galvannealed steel. Hot-dipped galvanized steel is annealed to create a diffusion alloyed layer.

Parameters

Grimm-Style DC 4 mm lamp
20 mA, 700 V
60 s at 100/s
180 s at 10/s

**For questions on this analysis e-mail us at:
info@leco.com**

For a complete listing of snapshots and performance notes visit us on the web at www.leco.com in the Spectroscopy section of the Applications Library.

