

Carbon/Nitrogen in Soil and Plant Tissue

LECO Corporation; Saint Joseph, Michigan USA

Instrument: TruMac[®] CN



Sample Preparation

Samples must be of uniform consistency to produce suitable results. Typically samples should be ground to a fineness of <0.5 mm.

Carbon and nitrogen results for soil and plant tissue materials are generally reported on a dry basis. The materials can either be dried prior to analysis or the moisture content determined and the values corrected. Please see the note addendum at the end of this document for details on drying these materials.

Accessories

528-203 Crucibles

Calibration Samples

502-092 EDTA, 502-642 Phenylalanine,
501-050 Nicotinic Acid

Analysis Parameters*

Furnace Temperature	1350°C
TE Cooler Temperature	5°C
Dehydration Time	0 seconds
Purge Cycles	3 seconds

Element Parameters

	Carbon	Nitrogen
Baseline Delay Time	0 seconds	6 seconds
Minimum Analysis Time	15 seconds	35 seconds
Comparator Level**	100.00	100.00
Endline Time	2 seconds	2 seconds
Conversion Factor	1.00	1.00
Significant Digits	5	5
TC Baseline Time	—	10 seconds
IR Analysis Stabilize Comparator	0	—
IR Baseline Time	1 second	—

Burn Profile

Burn Cycle	Lance Flow	Purge Flow	Time (seconds)
1	Off	On	5 seconds
2	On	On	5 seconds
3	On	Off	END

Ballast Parameters

Equilibrate Time	30 seconds
Not Filled Timeout	300 seconds

Aliquot Loop

Equilibrate Pressure Time	4 seconds
High Precision	Yes
High Speed	No

Procedure

1. Prepare instrument for operation as outlined in the operator's instruction manual.
2. Condition the system by analyzing 3 to 5 blanks (crucible is not required).
3. Determine blank.
 - a. Enter 1.0000 g mass into Sample Login (F3) using Blank as the sample name.
 - b. Place a 528-203 Crucible to the appropriate position of the autoloader.
 - c. Repeat steps 3a and 3b a minimum of three times.
 - d. Initiate the analysis sequence (F5).
 - e. Set the blank following the procedure outlined in the operator's instruction manual.
4. Calibrate.
 - a. Weigh ~0.50 g of EDTA calibration sample into a 528-203 Crucible, enter mass and sample identification into Sample Login (F3).
 - b. Transfer crucible to the appropriate position of the autoloader.
 - c. Repeat steps 4a and 4b a minimum of three times.
 - d. Initiate the analysis sequence (F5).
 - e. Calibrate the instrument following the procedure outlined in the operator's instruction manual.

Note: Multi-point (fractional weight or multiple calibration samples) may be used to calibrate if desired. A TruMac can be calibrated using several replicates of a single mass range (nominal 0.50 g) of EDTA utilizing a single standard calibration. The calibration can be verified by analyzing different compounds such as phenylalanine (0.25 g), nicotinic acid (0.15 g) or other reference materials with a similar matrix to the samples being analyzed.
5. Analyze Samples.
 - a. Weigh ~0.50 g sample into a 528-203 Crucible; enter mass and sample identification into Sample Login (F3).
 - b. Transfer crucible to the appropriate position of the autoloader.
 - c. Repeat steps 5a and 5b for each sample to be analyzed.
 - d. Initiate the analysis sequence (F5).

*Refer to TruMac Operator's Instruction Manual for Method Parameter definitions.

**Comparator Level indicator has been removed from methods as of software version 1.20.

Notes

Sample Drying Instructions:

- Soils – Samples dried at 105°C for one hour prior to analysis.
- Plant Tissues – Samples dried at 85°C for two hours prior to analysis.
- If soot (carbon black) is noticed in the primary filter (steel wool filter), reduce the sample mass to prevent soot build-up in this filter. Soot can be produced when larger masses of some sample types are analyzed.

Typical Results

(Based on a single standard calibration with 0.5 g of 502-092 EDTA)

Sample	Mass g	% Carbon	% Nitrogen	Sample	Mass g	% Carbon	% Nitrogen
Soil 502-308*	0.5155	2.80	0.234	Tobacco 502-082**	0.5021	46.55	2.54
Lot: 1010	0.5052	2.81	0.234	Lot: 1015	0.5055	46.51	2.54
2.77 ±0.06% C	0.5142	2.79	0.233	46.3 ±0.40% C	0.5001	46.56	2.54
0.233 ±0.013% N	0.5160	2.81	0.235	2.56 ±0.04% N	0.5065	46.71	2.54
	0.5005	2.82	0.236		0.5084	46.67	2.54
	0.5008	2.80	0.234		0.5020	46.59	2.54
	0.4998	2.82	0.235		0.5007	46.53	2.53
	0.5039	2.80	0.233		0.5007	46.49	2.53
	0.5006	2.82	0.235		0.5010	46.53	2.53
	0.5011	2.80	0.233		0.5025	46.47	2.52
	X =	2.81	0.234		X =	46.56	2.53
	s =	0.01	0.001		s =	0.07	0.01
Soil 502-309*	0.4989	13.34	1.15	Alfalfa 502-273**	0.5076	45.03	3.70
Lot: 1010	0.5010	13.35	1.14	Lot: 1021	0.5011	44.95	3.69
13.23 ±0.35% C	0.5051	13.21	1.13	45.0 ±0.3% C	0.5021	44.98	3.69
1.13 ±0.03% N	0.5003	13.39	1.14	3.74 ±0.03% N	0.5074	45.04	3.69
	0.4978	13.41	1.15		0.5055	44.90	3.68
	0.5007	13.28	1.14		0.5152	44.92	3.68
	0.5002	13.32	1.14		0.5345	44.86	3.68
	0.4996	13.31	1.14		0.5031	44.85	3.68
	0.5006	13.26	1.14		0.5063	44.87	3.68
	0.5034	13.35	1.15		0.5040	44.83	3.67
	X =	13.32	1.14		X =	44.92	3.68
	s =	0.06	0.01		s =	0.07	0.01
Orchard Leaves	0.5025	50.99	2.39				
502-055**	0.5011	50.89	2.38				
Lot: 1031	0.5032	50.78	2.36				
50.7 ±0.3% C	0.5048	51.00	2.37				
2.39 ±0.06% N	0.5013	50.83	2.36				
	0.5037	51.02	2.37				
	0.5018	50.88	2.36				
	0.5110	50.84	2.36				
	0.5019	50.83	2.36				
	0.5047	50.91	2.36				
	X =	50.90	2.37				
	s =	0.08	0.01				

*Samples were dried at 105°C for one hour prior to analysis.

**Samples were dried at 85°C for two hours prior to analysis.

TruMac is a registered trademark of LECO Corporation.

LECO Corporation | 3000 Lakeview Avenue | St. Joseph, MI 49085 | Phone: 800-292-6141 | Fax: 269-982-8977
 info@leco.com • www.leco.com | ISO-9001:2008 HQ-Q-994 | LECO is a registered trademark of LECO Corporation.