


Note: Part numbers and standards' values may change. Consult LECO for the latest information.



See reverse side for contact information

GENERAL COMBUSTION SUPPLIES

501-171-HAZ Anhydron 454 g (1 lb.) bottle
 502-190 Copper Oxide Wire 100 g bottle
 502-705 Copper Sticks 100 g ampule
 502-878 Copper Sticks 100 g bottle
 502-656 Copper Turnings 60 g bottle
 501-081 Glass Wool 454 g (1 lb.) box
 763-265 Glass Wool 227 g (0.5 lb.) package
 502-174-HAZ Lecosorb 500 g bottle
 502-644 Micro Furnace Reagent
 502-049 N-Catalyst 50 g bottle
 608-379 Quartz Wool Strips 10/pkg
 502-301 Tungsten Trioxide (125 g)

CAPSULES AND CRUCIBLES

614-961-110 Porous Crucibles for solids 10/pkg (use only with 614-962 lance)
 502-293 Silver Capsules for Micro Analysis 500/bottle
 502-206 Silver Capsules for Micro Analysis 1000/bottle
 502-227 Tin Capsules for Micro Analysis 500/bottle
 502-489 Reticulated Micro Crucibles 10/pkg (for CHNS only)
 625-511-154 Porous Reticulated Crucibles (small) 5/pkg

OPERATION ACCESSORIES

502-005 Brush	604-398 Sample Cup Holder, Large
601-981 Wire Brush	604-373 Sample Cup Holder, Small
616-152 Crucible Extractor Tool	778-405 10-Sample Storage Block
775-306 Secondary Filter, 10 micron	501-614 Spatula
775-307 Filter Tool, Thumbscrew	502-007 Tube Cleaner
502-023 Funnel Quick Disc	502-028 Tweezers
778-321 Lance Extraction Tool	604-378 Tweezers, Curved
502-213 Leak Detection Solution	501-241 Vacuum Grease
603-350 Tool O-ring Installation	619-380 Vacuum Cleaner Attachment Assembly
762-515 Tool O-ring Removal	

COMBUSTION TUBES AND GLASSWARE

619-701 Combustion Tube, Micro
 619-425 Micro In-Line Anhydron Tube
 774-638 Reagent Tube, requires 772-520 O-rings for sealing
 775-601 Reagent Tube, requires 772-520 O-rings for sealing
 619-154 Reduction Tube, requires 601-442 O-rings for sealing

SET-UP STANDARDS**EDTA (50 grams)**

Part No.	% C	% H	% N	% S
502-092	40.80–41.20	5.46–5.54	9.52–9.60	—

SUCROSE (50 grams)

Part No.	% C	% H	% N	% O	% S
501-441	42.10	6.48	—	—	—

POTASSIUM BIPHTHALATE (20 grams)

Part No.	% C	% H	% N	% O	% S
501-048	47.05	2.47	—	—	—

ACETANILIDE (10 grams)

Part No.	% C	% H	% N	% O	% S
501-053	71.09	6.71	10.36	11.84	—

BENZOIC ACID (10 grams)

Part No.	% C	% H	% N	% O	% S
502-184	68.85	4.92	—	26.23	—

BBOT (25 grams)

Part No.	% C	% H	% N	% O	% S
502-654	72.5	6.0	6.5	—	7.4

DURENE (10 grams)

Part No.	% C	% H	% N	% O	% S
502-634	89.49	10.51	—	—	—

STEARIC ACID (5 grams)

Part No.	% C	% H	% N	% O	% S
502-203	76.00	12.76	—	—	—

DINITROBENZOIC ACID (5 grams)

Part No.	% C	% H	% N	% O	% S
502-204	40.39	1.92	11.54	46.15	—

CAFFEINE (5 grams)

Part No.	% C	% H	% N	% O	% S
502-205	49.48	5.15	28.87	16.50	—

CYSTINE (5 grams)

Part No.	% C	% H	% N	% O	% S
502-207	29.99	5.03	11.65	26.65	26.68

PHENYLALANINE (5 g)

Part No.	% C	% H	% N	% O	% S
502-641	65.4	6.7	8.5	—	—

SULFAMETHAZINE (5 grams)

Part No.	% C	% H	% N	% O	% S
502-209	51.78	5.07	20.13	11.5	11.52

*Values vary—contact LECO for current levels.

NOTE: Part numbers and standards' values may change.

Consult LECO for latest information.